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Grating Structured Surface Plasmon Resonance Chip and Plasmon Waveguide Resonance/Surface Plasmon Resonance Combined Chip

PI: Prof. Ling Chao

Department of Chemical Engineering,


National Taiwan U.

Experience:

Website: http://web.che.ntu.edu.tw/biomembrane/

Market Needs:

Recently, lots of researchers use surface plasmon resonance and its related technique to

study in the fields of biology and chemistry. The developed chip in this study provides more information than conventional chip did, and we demonstrate the new application for surface plasmon resonance. Newly developed chip can be applied directly in current commercial equipment without the change or adjustment in parts.

Our Technology:

In this study, we created PWR/SPR combined chip that it can retain both surface plasmon resonance (SPR) and plasmon-waveguide resonance (PWR). This chip has sub-micron sized grating pore array structures. In the pores, resonance film is coated at the bottom and the change of refractive index in the pore region close to the surface of resonance film is sensed by the SPR. In the place without etching, PWR occurs in resonance film with waveguide layer. The geometry of the PWR/SPR combined chip allowed us to use the SPR to detect the refractive index change in the pore region, which is correlated to the target species concentration inside the pore, and the PWR to simultaneously monitor the change of refractive index at the top waveguide layer surface. This chip provides two resonance peaks that can react to change in different area, respectively.

Strength:

- 1. Provide two resonace peaks that react to changes happen in different area.
- 2. It can be used in current commercial machine without adjusting the parts.
- 3. Sense the physical properties so that it can be cooperated with current chemical detection platform.

Competing Products:

- 1. Conventional SPR chip. However, it has one peak and limited sensing area.
- 2. PWR chip. However, it has one peak and limited sensing area.
- 3. LSPR chip. However, it has one peak and limited sensing area.

This information herein is intended for potential license of NTU technology only. Other usage of all or portion of this information in whatever form or means is strictly prohibited. Kindly contact us and we will help to achieve your goal the best we can.

Intellectual Properties:

Our group have the techniques: chip fabrication and design, optical system design, cell culture and membrane study, operation and design of microfluidics system and simulation. These techniques are related application to support the further development of surface palsmon resonance.