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Marketing Abstract of NTU's Invention Disclosure

NTU's docket no: 09T-100219

CIAC contact : Lou Wei-Chen

Tel : 02-33669948

e-mail : weichenlou@ntu.edu.tw

Title	<p>There are three products.</p> <ol style="list-style-type: none"> 1. CD8α cDNA plasmid. CD8α was the cell surface marker of grouper (<i>Epinephelus coioides</i>) cytotoxic T lymphocyte. 2. CD8α recombinant protein 3. Rabbit antiserum against CD8α.
Inventor (s)	<p>Professor Yen-Ling Song, Institute of Zoology and Department of Life Science, College of Life Science Yun-Tsan Chang, Institute of Zoology</p>
Brief Description	<ol style="list-style-type: none"> 1. CD8α cDNA was cloned and characterized from grouper (<i>Epinephelus coioides</i>) cytotoxic T lymphocyte membrane protein 2. Recombinant protein of grouper CD8α was biosynthesized using the insect cell-baculovirus expression system. 3. Rabbit antiserum was prepared against CD8α recombinant protein at 8-9 weeks post immunization. 4. The antiserum recognized specifically the subtype of peripheral lymphocytes as confirmed by immunocytochemical stain. 5. The ELISA titer can reach to 10,000 x and the ICC titer to 5,000 x.
Fields of Application	<ol style="list-style-type: none"> 1. The plasmid containing CD8α cDNA can express the CD8α recombinant protein. 2. The CD8α recombinant protein can be used for preparation of the CD8α antibody. The CD8α recombinant protein can also serve as a competitor for CTL binding to MHC class I molecule in the immunoassay. 3. The CD8α antiserum/antibody can be applied in ELISA 、Western blot 、immunocytochemical (ICC) and immunohistochemical (IHC) stains.
Advantages	<ol style="list-style-type: none"> 1. There are few antibodies available for fish commercially, particularly for fish leucocyte or lymphocyte antibodies. 2. This product will be an important tool for studying the cellular immune responses in fish and useful for evaluating the fish vaccine efficacy. 3. This is the sole product of grouper cytotoxic T lymphocyte CD8α antibody at present.
Market Potential	<ol style="list-style-type: none"> 1. Grouper (<i>Epinephelus coioides</i>) is an important fish species for Asian aquaculture. 2. CD8α antibody is an important tool for studying the cellular immune responses in fish.
IP Right(s)	<p>Know-how</p>

Picture


Fig. 1. (A) Reducing SDS-PAGE of recombinant protein stained with Coomassie blue. The arrow head indicates the target product of MW 24 kDa. (B) Western blot of the target product was recognized by anti-6xHis tag antibodies (arrow head). MM: MW marker; E1: Eluted product 1; E2: Eluted product 2.


Fig. 2. ELISA titers of rabbit antisera against recombinant grouper CD8 α in the weeks post immunization. Rabbit antiserum was 2-fold serially diluted.


Fig. 3. Confocal microscopy of grouper peripheral blood leucocytes stained with rabbit antiserum against recombinant CD8 α (W8, 500x dilution). Rabbit antiserum against grouper CD8 α recombinant protein served as the first antibody (W8, 500x). FITC-conjugated goat anti-rabbit Ig served as the secondary antibody. Characteristic round- or ovoid-shaped cells with central monomorphonucleus, high nuclear to cytoplasmic ratio and ring-form cytoplasm indicate the subset lymphocytes. Bar =25 μ m.

Picture


Fig. 4. Specific cytotoxicity of CD8 α^+ PBL effectors. Harvested from fish at one week post-NNV sensitization, CD8 α^+ PBL effectors enhanced the cytotoxicity statistically significantly against NNV-infected targets in comparison with RSIV-infected targets. The specific cytotoxicity of CD8 α^+ PBLs revealed that the lymphocyte subset recognized by CD8 α antibodies are the CTLs. Data represent mean \pm S.D. (n=3).